


MS Spectrum



Sample Information :	Interface :	ESI	Perod Bias :	+4.5kv
Order ID: :	Nebulizing Gas Flow :	1.5L/min	Detector :	-0.2kv
Dissolution method :	CDL Temp :	250C	T.Flow :	0.2ml/min
Date Acquired :	CDL Volt :	0v	B.conc :	30%H ₂ O/70%MeOH
Injection Volume :	Block Temp :	200		

Lot No.	TB400451- TB400784
Theoretical	4963.5
Observed	4962.9